

# **Bed Rail Guide**

Rails options and configurations compatible with our different bed styles

## **Full Electric Bed Frame 36" (Invacare/Drive)**

- 0 Rails
- 1 Half Length Side Rails (Pair) Cannot use second pair on the same bed for safety reasons. See Bari beds if needing that 4 rail arrangement.
- 1 Full Length Side Rails (Pair)
- 1 Assist Side Rails (Pair)
- 1 Halo Rails (Pair) Check availability

## **Full Electric Low Bed 36" (Invacare/Drive)**

- 0 Rails
- 1 Assist Side Rails (Pair) Cannot use half or full rails on this model for safety reasons.
- 1 Halo Rails (Pair) Check availability

#### FloorBed Frame 36"or 42"

- 0 Rails
- 1 FloorBed Assist Side Rail (Single) One side
- 2 FloorBed Assist Side Rail (Pair) Both sides

### **Bariatric Bed Frame 42" or 48"**

- 0 Rails
- 1 Bariatric Side Rails (Pair)
- 1 Bariatric Side Rails (4)

## **Bariatric Bed Frame (MedaCure Only)**

- 0 Rails
- 1 Bariatric Side Rails (Pair)
- 1 Bariatric Side Rails (4)
- 1 Halo Rails (Pair) Specifically ok on this model. Not compatible with other bari beds.