

# Hemo-Force DVT Prophylaxis Pumps

Medline Hemo-Force and Hemo-Force II pumps offer both intermittent and sequential options to address clinician judgment, patient needs and facility protocols. They are suitable for DVT reduction in acute care environments, in the pre-, intra-, and post-surgery phases, as mono therapy or when combined with anti-coagulation therapy.



**NEW!**

## Hemo-Force II

**25% Smaller and 100% Auto-Detection**

- » Allows a combination of sleeves to be used simultaneously: foot and calf, foot and thigh, calf and thigh
- » Auto detects quantity and type of garment without any additional user input
- » Compact, durable design with integrated handle for easy transport

### Specifications

Screen	LCD with real-time display
Timer	Minutes, hours
Alarm	Audible, visual, reference code and mute button
Tubing	Anti-kink, easy to clean, loss protection
Hooks	Spring-loaded, fold-away, rubberized and metal reinforced
Troubleshooting Guide	Tri-lingual
Therapy Types Available	Intermittent and sequential (separate pumps), single or dual leg therapy
Sleeve Types Compatible	Foot/calf/thigh; any combination
Inflation Time	Intermittent: 12 seconds Sequential: 12 seconds (4 seconds per bladder)
Deflation Time	48 seconds
Battery	Lithium ion (optional, sold separately)
Power Cord	2-Prong, 15-ft. (Class II, double insulated)
Weight	5-lbs.

### Hemo-Force II DVT Pumps

Item No.	Description	Pkg.
MDS600INT2	Hemo-Force II Intermittent	1/ea
MDS600SQ2	Hemo-Force II Sequential	1/ea
MDS600BAT	DVT Pump Battery	1/ea